



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Information

DEAC-dUTP, lyophilized powder

Catalog Number 40059

Unit Size 25 nmole

Molecular Weight 762

Color and Form Light yellow solid

Absorption/Emission 426/480 nm

Storage and Handling

Store DEAC-dUTP at -20°C or -70°C, protected from light. When stored as recommended, product is stable for at least 6 months from date of receipt.

Product Description

Fluorescent conjugates of dUTP can be used for detection of apoptotic cells by direct fluorescence TUNEL labeling of DNA strand breaks in cells¹, and can be used in place of dTTP in standard DNA labeling and synthesis protocols to generate fluorescent dsDNA and oligonucleotide probes.

Note: The product is not compatible with PCR applications.

References

1. Gold et al. (1994). *Lab Invest.* 71 (2):219-25.
2. Slupphaug et al. (1993). *Anal Biochem.* 211 (1):164-9.
3. Hogrefe et al. (2002). *PNAS* 99 (2): 596-601.

Please visit our website at www.biotium.com to view our full selection of CF™ dye bioconjugates, including antibodies, antibody labeling kits, phalloidin, Annexin V and α -bungarotoxin, as well as classic fluorescent nucleic acid dyes and antifade fluorescence mounting media.