

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

Product Information



CRTH2/DP, Receptor (N-Term) Blocking Peptide

Item No. 10004884

Prostaglandin D2 (PGD2) elicits its biological function through interaction with two distinct G protein-coupled receptors, DP1 and CRTH2/DP2. CRTH2 primarily couples to G1 or G1 subunits to mobilize Ca2+, induce cell migration, and up-regulate adhesion molecules.¹ CRTH2 mRNA has been detected in various tissues including liver, lung, kidney, brain, heart, thymus, and spleen and in various cell lineages including both hematopoietic and non-hematopoietic cell lines.² Human CRTH2 is 395 amino acids in length with an estimated molecular weight of 43 kDa. Cayman's CRTH2/ DP2 receptor (N-Term) polyclonal antibody can be used for western blot and immunohistochemical analysis of CRTH2 on samples of human, bovine, murine, and rat origin. It detects both unglycosylated and glycosylated protein at sizes ranging from 35-40 kDa to 50-70 kDa, as reported by Nagata et al., in 1999.³

Laboratory Procedures

This vial contains 200 μg peptide in 200 μl TBS, pH 7.4, containing 0.1% BSA and 0.02% sodium azide. The CRTH2/DP₂ receptor (N-Term) blocking peptide (human CRTH2 amino acids 2-21) can be used in conjunction with Cayman's CRTH2/DP2 Receptor (N-Term) Polyclonal Antibody (Item No. 10004886) to block protein-antibody complex formation during immunochemical analysis of CRTH2.

Store this peptide solution at -20°C. It will be stable for at least two years. To block antibody/protein complex formation, the following procedure is recommended:

- 1. Mix the CRTH2/DP2 Receptor (N-Term) Polyclonal Antibody (Item No. 10004886) and blocking peptide together in a 1:1 (v/v) ratio in a microfuge tube. For example, mix 20 µl of antibody and 20 µl of peptide.*
- 2. Incubate for one hour at room temperature with occasional mixing prior to further dilution and application of the mixture to the immunoblot.

3. Dilute the mixture to the final working antibody concentration and apply to the slide or membrane as usual.

*This is a recommended mixture. The minimum amount of peptide needed for complete blocking has not been precisely determined and may vary depending on the sample being analyzed. The amount of peptide required may need to be increased if sufficient blocking does not occur.

References

- 1. Nagata, K. and Hirai, H. The second PGD₂ receptor CRTH2: Structure, properties, and functions in leukocytes. Prostaglandins Leukot. Essent. Fatty Acids 69, 169-177 (2003).
- Abe, H., Takeshita, T., Nagata, K., et al. Molecular cloning, chromosome mapping and characterization of the mouse 2. CRTH2 gene, a putative member of the leukocyte chemoattractant receptor family. Gene 227, 71-77 (1999).
- 3. Nagata, K., Tanaka, K., Ogawa, K., et al. Selective expression of a novel surface molecule by human Th2 cells in vivo. J. Immunol. 162, 1278-1286 (1999).

Related Products

For a list of related products please visit: www.caymanchem.com/catalog/10004884

Cayman Chemical

Mailing address 1180 E. Ellsworth Road Ann Arbor, MI 48108 USA

Phone (800) 364-9897 (734) 971-3335

Fax (734) 971-3640

E-Mail custserv@caymanchem.com

Web www.caymanchem.com

WARNING: THIS PRODUCT IS FOR LABORATORY RESEARCH ONLY: NOT FOR ADMINISTRATION TO HUMANS. NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

MATERIAL SAFETY DATA

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user <u>must</u> review the <u>complete</u> Material Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular prose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman warrants only to the original customer that the material will meet our specifications

purpose, suitability and merchantability, which extends beyond the description of the chemicals nervor. Cayman warrants only to the original customer that the material will <u>meet our specifications</u> at the time of delivery. Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, indirect, incidental or consequential damages, even if Cayman is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees. Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a <u>refund</u> of the purchase price, or at Cayman's option, the <u>replacement</u>, at no cost to Buyer, of all material that does not meet our specifications. Said refund as repleciment is conditioned on Buyer a pointe to common within thirty (20) days after a return of the purchase of Buyer to noise point pointer of the purchase of Buyer to noise point pointer of the purchase of the purchase of the purchase of the purchase of Buyer to noise point point point of the purchase of Buyer to point point point point of the purchase of Buyer to point point point of the purchase of Buyer to point po

es not meet our specifications. Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within try (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material. For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog. Copyright Cayman Chemical Company, 07/18/2012