



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

PRODUCT INFORMATION



Cell-Based Assay IBMX Solution (1,000X)

Item No. 10008978

Overview and Properties

Contents: This vial contains 500 µl of 1,000X IBMX stock solution.
Storage: -20°C
Stability: ≥3 years

Procedures

To prepare 10 ml of induction culture medium, add 10 µl of this IBMX (1,000X) solution, along with 10 µl of Insulin Solution (1,000X) (Item No. 10008979) and 10 µl of Adipogenesis Assay Dexamethasone Solution (1,000X) (Item No. 10008980) to DMEM containing 10% FBS. To prepare 10 ml of 2X IBMX Solution for STEP Plate Assay Kits, add 10 µl of this stock solution to the culture medium used for the cells. To avoid frequent freeze and thaw of the stock solution, we recommend that you make small aliquots for long-term storage.

CAUTIONS

For laboratory research use only. This reagent is not for human or veterinary use. This product has been tested and formulated to work exclusively with the Adipogenesis Assay Kit and STEP Plate Assay Kits (cAMP Method) from Cayman (Item Nos. 10006908, 600380, 600390, 600400, 600410, 600420, and 600430). This product may not perform as described if used with other assay reagents or protocols.

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 10/25/2021

CAYMAN CHEMICAL
1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA
PHONE: [800] 364-9897
[734] 971-3335
FAX: [734] 971-3640
CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM