

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# **Product Information**



## Nuclear Extraction Buffer (2X)

Catalog No. 10009306

Contents: This vial contains 3 ml of 2X buffer.

Stability: ≥1 year at 4°C

Application: For use with the nuclear extraction kit.

## **Laboratory Procedures**

To prepare the buffer for use in the Nuclear Extraction Kit (Catalog No. 10009277), dilute reagents as outlined below. Some of the phosphatase and protease inhibitors lose activity shortly after dilution, therefore any remaining 2X buffer should be discarded.

Reagent	60 mm plate*	100 mm plate**	150 mm plate***
Nuclear Extraction Buffer (2X)	12 μl	25 μl	50 μl
Nuclear Extraction Protease Inhibitors (100X)	0.25 μl	0.5 μl	1.0 µl
Phosphatase Inhibitor Cocktail (50X)	0.5 μl	1.0 µl	2.0 µl
Nuclear Extraction DTT (10 mM)	2.5 µl	5.0 μl	10 μl
Distilled water	9.25 µl	18.5 µl	37 μl
Total Volume	25 μl	50 μl	100 µl

<sup>\*60</sup> mm plate is  $-3.5 \times 10^6$  cells, \*\*100 mm plate is  $-7 \times 10^6$  cells, \*\*\*150 mm plate is  $-1.5 \times 10^7$  cells

For laboratory research use only. This reagent is not for human or veterinary use. This product has been tested and formulated to work exclusively with the Nuclear Extraction Kit from Cayman (Catalog No. 10009277). This product may not perform as described if used with other assay reagents or protocols.

## **Related Products**

Nonidet P-40 Assay Reagent (10%) - Cat. No. 600009 • Nuclear Extraction Kit - Cat. No. 10009277 • Nuclear Extraction Hypotonic Buffer (10X) - Cat. No. 10009301 • Nuclear Extraction Dithiothreitol (1M) - Cat. No. 10009302 • Nuclear Extraction Protease Inhibitor Cocktail (100X) - Cat. No. 10009303 • Nuclear Extraction PBS (10X) - Cat. No. 10009304 • Nuclear Extraction Phosphatase Inhibitors (50X) - Cat. No. 10009305

WARNING: This product is for laboratory research only: not for administration to humans. Not for human or veterinary DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user must review the complete Material Safety Data

WARRANTY AND LIMITATION OF REMEDY Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular pose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman warrants only to the original customer that the material will meet our specifications

at the time of delivery.

Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have any obligation or liability, whether in tort (including negligence) or in contract, for

direct, indirect, incidental or consequential damages, even if Cayman is informed about their possible existence.

This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees.

Buyer's exclusive remedy and Cayman's sole liability hereunder shall be limited to a refund of the purchase price, or at Cayman's option, the replacement, at no cost to Buyer, of all material that

does not meet our specifications.

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.

Copyright Cayman Chemical Company, 01/21/2010

## Cayman Chemical

## **Mailing address**

1180 E. Ellsworth Road Ann Arbor, MI 48108 USA

## **Phone**

(800) 364-9897 (734) 971-3335

(734) 971-3640

custserv@caymanchem.com

www.caymanchem.com