



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Product Information



## SIRT1 FRET Developer

Item No. 10011085

**Contents:** The vial contains 50 µl of a lysyl endopeptidase solution.

**Stability:** ≥1 year at -80°C

### Laboratory Procedures

Thaw the developer on ice. Dilute 15 µl of developer with 735 µl of diluted SIRT1 FRET Assay Buffer (10X) (Item No. 10011080) prior to performing the SIRT1 FRET-Based Screening Assay Kit (Item No. 10010991). The diluted developer is stable for four hours at 4°C.

### CAUTIONS

For laboratory research use only. This reagent is not for human or veterinary use. This product has been tested and formulated to work exclusively with the SIRT1 FRET-Based Screening Assay Kit from Cayman (Item No. 10010991). This product may not perform as described if used with other assay reagents or protocols.

### Related Product

SIRT1 FRET-Based Screening Assay Kit - Item No. 10010991

## Cayman Chemical

### Mailing address

1180 E. Ellsworth Road  
Ann Arbor, MI  
48108 USA

### Phone

(800) 364-9897  
(734) 971-3335

### Fax

(734) 971-3640

### E-Mail

custserv@caymanchem.com

### Web

www.caymanchem.com

**WARNING: THIS PRODUCT IS FOR LABORATORY RESEARCH ONLY: NOT FOR ADMINISTRATION TO HUMANS. NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

#### MATERIAL SAFETY DATA

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user must review the complete Material Safety Data Sheet, which has been sent *via* email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman **warrants only** to the original customer that the material will meet our specifications at the time of delivery.

Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, indirect, incidental or consequential damages, even if Cayman is informed about their possible existence.

This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees.

Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a refund of the purchase price, or at Cayman's option, the replacement, at no cost to Buyer, of all material that does not meet our specifications.

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

**For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.**