



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Product Information



## Sampling tubes with Indomethacin

Cat No: D31008 - 35 tubes

### General Data

- Shipping:** room temp.
- Stability:** Store at RT; shelf life 2 years maximum after production



- Application(s):** Indomethacin sampling tubes are ready to use. Just collect your samples in classical blood collection tubes, and transfer 1 mL of sample per tube.

These tubes are not intended to be used to directly collect blood.

### Product Overview

Indomethacin is a cyclooxygenase (COX) inhibitor that prevents ex-vivo Thromboxane B2 formation and other prostaglandins that can not be prevented using anticoagulant alone. The samples have to be transferred in Indomethacin sampling tubes as soon as possible after blood collection. For measuring some prostaglandins, it may be useful to block the ex-vivo formation of prostaglandin. The tubes are graduated, with a silicone washer sealed and external threads. They can be stored at temperatures as low as -80 degrees. Thus, indomethacin sampling tubes may be useful for any COX metabolites studies.

### Scientific Literature

Ferreira S, Moncada S, Vane J (Jun 23 1971). "Indomethacin and aspirin abolish prostaglandin release from the spleen". Nat New Biol 231 (25): 237-9

FP/08/17

#### For research laboratory use only - Not for human diagnostic use.

Buyers agree to purchase the material subject to Purchasing Terms that can be found on our website. Seek appropriate training to safely handle this product under normal conditions. Use the recommended personal protective equipment to prevent chemical exposures.

Bertin Pharma does not make any other warranty or representation whatsoever whether expressed or implied, with respect to these products. In no event will Bertin Pharma be liable for incidental, consequential or punitive damages.

#### Contact Bertin Pharma

10 bis avenue Ampère  
F-78180 Montigny-le-Bretonneux  
Tel.: +33 (0)139 306 036  
bioreagent@bertinpharma.com  
**bioreagent.bertinpharma.com**