

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



7-keto Cholesterol

Item No. 16339

CAS Registry No.: 566-28-9

Formal Name: 3β-hydroxy-cholest-5-en-7-one Δ^5 -Cholesterol-3 β -ol-7-one, Synonyms:


7-oxo Cholesterol, SC-4722

MF: $C_{27}H_{44}O_2$ FW: 400.6 **Purity:** ≥98% UV/Vis.:

 λ_{max} : 238 nm A crystalline solid Supplied as:

-20°C Storage: Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



Laboratory Procedures

7-keto Cholesterol is supplied as a crystalline solid. A stock solution may be made by dissolving the 7-keto cholesterol in the solvent of choice. 7-keto Cholesterol is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide, which should be purged with an inert gas. The solubility of 7-keto cholesterol in these solvents is approximately 20, 0.1, and 2 mg/ml, respectively.

7-keto Cholesterol is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, 7-keto Cholesterol should first be dissolved in ethanol and then diluted with the aqueous buffer of choice. 7-keto Cholesterol has a solubility of approximately 0.3 mg/ml in a 1:2 solution of ethanol:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

7-keto Cholesterol is a bioactive sterol and a major oxysterol component of oxidized LDL.^{1,2} It is produced by oxidation of cholesterol via ethanol-mediated lipid peroxidation or photodamage as well as oxidation of 7-dehydro cholesterol (Item No. 14612) by the cytochrome P450 (CYP) isoform CYP7A1.3-5 7-keto Cholesterol inhibits CYP7A1 (IC $_{50}$ = ~1 μ M). It induces activation and chemotaxis of retinal microglia as well as polarization to a pro-inflammatory state via NLRP3 inflammasome activation in vitro.⁶ Intraocular implantation of 7-keto cholesterol coated wafers increases ocular levels of VEGF, IL-1β, and GRO/KC, macrophage infiltration, and neovascularization in rat eye. Levels of 7-keto cholesterol in lipid deposits are increased in a variety of chronic diseases, including atherosclerosis, Alzheimer's disease, and age-related macular degeneration.

References

- 1. Shentu, T.P., Titushkin, I., Singh, D.K., et al. Am J Physiol Cell Physiol. 299(2), C218-C229 (2010).
- 2. Brown, A.J., Leong, S.I., Dean, R.T., et al. J. Lipid Res. 38(9), 1730-1745 (1997).
- 3. Rodriguez, I.R. and Fliesler, S.J. Photochem. Photobiol. 85(5), 1116-1125 (2009).
- 4. Shinkyo, R., Xu, L., Tallman, K.A., et al. J. Biol. Chem. 286(38), 33021-33028 (2011).
- Mitić, T., Shave, S., Semjonous, N., et al. Biochem. Pharmacol. 86(1), 146-153 (2013).
- 6. Indaram, M., Ma, W., Zhao, L., et al. Sci. Rep. 5:9144, (2015).
- 7. Amaral, J., Lee, J.W., Chou, J., et al. PLoS One 8(2), e56099 (2013).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after ndling. Before use, the user <u>must</u> review the <u>complete</u> Safety Data Sheet, which has been sent *via* email to yo

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information

Copyright Cayman Chemical Company, 08/30/2018

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA **PHONE:** [800] 364-9897

[734] 971-3335 **FAX:** [734] 971-3640 CUSTSERV@CAYMANCHEM.COM

WWW.**CAYMANCHEM**.COM