

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Glycohyocholic Acid

Item No. 22670

CAS Registry No.: 32747-08-3

Formal Name: N-[$(3\alpha,5\beta,6\alpha,7\alpha)$ -3,6,7-trihydroxy-

24-oxocholan-24-yl]-glycine

Synonyms: GHCA, Glycine Hyocholate,

Glycohyocholate

MF: $C_{26}H_{43}NO_{6}$ FW: 465.6 **Purity:** ≥98%

UV/Vis.: λ_{max} : 202 nm Supplied as: A crystalline solid

-20°C Storage: Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

ÓН

Laboratory Procedures

Glycohyocholic acid (GHCA) is supplied as a crystalline solid. A stock solution may be made by dissolving the GHCA in the solvent of choice. GHCA is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF), which should be purged with an inert gas. The solubility of GHCA in ethanol is approximately 20 mg/ml and 30 mg/ml in DMSO and DMF.

GHCA is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, GHCA should first be dissolved in DMSO and then diluted with the aqueous buffer of choice. GHCA has a solubility of approximately 0.5 mg/ml in a 1:1 solution of DMSO:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

GHCA is a mammalian secondary bile acid and glycine conjugate of hyocholic acid (Item No. 20293).¹ GHCA is upregulated 6.6-fold in the serum of patients with cirrhosis induced by hepatitis C virus (HCV) compared to healthy controls.² Plasma levels of GHCA increase in patients who are no longer diabetic following gastric bypass surgery.3

References

- 1. Guo, L., Milburn, M.V., Ryals, J.A., et al. Plasma metabolomic profiles enhance precision medicine for volunteers of normal health. Proc. Nat. Acad. Sci. USA 112(35), E4901-E4910 (2015).
- 2. Fitian, A.I., Nelson, D.R., Liu, C., et al. Integrated metabolomic profiling of hepatocellular carcinoma in hepatitis C cirrhosis through GC/MS and UPLC/MS-MS. Liver Int. 34(9), 1428-1444 (2014).
- 3. Yu, H., Ni, Y., Bao, Y., et al. Chenodeoxycholic acid as a potential prognostic marker for roux-en-Y gastric bypass in Chinese obese patients. J. Clin. Endocrinol. Metab. 100(11), 4222-4230 (2015).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

al should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution

WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 01/23/2018

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA **PHONE:** [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM