



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 


# PRODUCT INFORMATION



## Xanthine amine congener

Item No. 23077

**CAS Registry No.:** 96865-92-8  
**Formal Name:** N-(2-aminoethyl)-2-[4-(2,3,6,7-tetrahydro-2,6-dioxo-1,3-dipropyl-1H-purin-8-yl)phenoxy]-acetamide  
**Synonym:** XAC  
**MF:** C<sub>21</sub>H<sub>28</sub>N<sub>6</sub>O<sub>4</sub>  
**FW:** 428.5  
**Purity:** ≥98%  
**Supplied as:** A solid  
**Storage:** -20°C  
**Stability:** ≥2 years



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

### Laboratory Procedures

Xanthine amine congener (XAC) is supplied as a solid. A stock solution may be made by dissolving the XAC in the solvent of choice. XAC is soluble in the organic solvent DMSO, which should be purged with an inert gas.

### Description

XAC is a non-selective adenosine receptor antagonist with  $K_b$  values of 83, 25, and 15 nM for rat PC12 cells, human platelets that endogenously express adenosine A<sub>2A</sub> receptors, and rat fat cells that endogenously express A<sub>1</sub> receptors, respectively.<sup>1</sup> It reduces adenylate cyclase activity stimulated by 5'-N-ethylcarboxamidoadenosine (NECA; Item No. 21420) in human platelets ( $EC_{50}$ s = 1 and 0.31 μM in the presence and absence of XAC, respectively). XAC also reverses R-PIA-induced inhibition of isoproterenol-stimulated adenylate cyclase activity in rat fat cell membranes ( $IC_{50}$ s = 146 and 26 nM in the presence and absence of XAC, respectively). *In vivo*, XAC reverses reductions in heart rate and blood pressure induced by 2-chloro-3-deazaadenosine (2-CADO; Item No. 9000787), R-PIA, and NECA in rats when administered at doses ranging from 0.1 to 1.0 mg/kg.<sup>2</sup> XAC (39.8 mg/kg) also induces convulsions in mice.<sup>3</sup>

### References

1. Ukena, D., Daly, J.W., Kirk, K.L., *et al.* Functionalized congeners of 1,3-dipropyl-8-phenylxanthine: Potent antagonists for adenosine receptors that modulate membrane adenylate cyclase in pheochromocytoma cells, platelets and fat cells. *Life Sci.* **38(9)**, 797-807 (1986).
2. Evoniuk, G., Jacobson, K.A., Shamim, M.T., *et al.* A1- and A2-selective adenosine antagonists: *In vivo* characterization of cardiovascular effects. *J. Pharmacol. Exp. Ther.* **242(3)**, 882-887 (1987).
3. Morgan, P.F., Deckert, J., Jacobson, K.A., *et al.* Potent convulsant actions of the adenosine receptor antagonist, xanthine amine congener (XAC). *Life Sci.* **45(8)**, 719-728 (1989).

#### WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

#### SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 11/18/2020

#### CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD  
ANN ARBOR, MI 48108 · USA

**PHONE:** [800] 364-9897

[734] 971-3335

**FAX:** [734] 971-3640

CUSTSERV@CAYMANCHEM.COM  
WWW.CAYMANCHEM.COM