

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



20-HETE-d₆ MaxSpec[®] Standard *Item No. 24923*

Formal Name:	20-hydroxy-5Z,8Z,11Z,14Z- eicosatetraenoic-16,16,17,17,18,18-d ₆ acid соон
Synonym:	20-hydroxy Arachidonic Acid-d ₆ $/$ $/$ $/$ $/$ $/$ D D D
MF:	$C_{20}H_{26}D_6O_3$
FW:	326.5 ОН
Purity:	≥95% D D
Supplied as:	A solution in acetonitrile; in a deactivated glass ampule
Concentration:	10 μg/ml (nominal); see certificate of analysis for verified concentration
Storage:	-20°C
Stability:	≥3 years; Stability testing is ongoing to ensure concentration accuracy. The certificate of analysis and
	product expiry date will be updated upon completion of testing.
Special Conditions: Store upright and unopened at -20°C. Warm to room temperature prior to opening.	
	Light sensitive.

Description

20-HETE-d₄ is intended for use as an internal standard for the quantification of 20-HETE (Item No. 90030) by GC- or LC-MS. 20-HETE is a cytochrome P450 (CYP450) metabolite postulated to play an autacoid role in the renal and cerebral vasculature.¹ In rat cerebral microvessels, 20-HETE is a vasoconstrictor that mediates pressure-induced autoregulatory vasoconstriction.² 20-HETE is excreted mainly as the glucuronide conjugate. The concentration of free 20-HETE (20-40 pg/ml in human urine) is about 10-fold lower than the corresponding concentration of the 20-glucuronide.³

20-HETE-d, MaxSpec[®] standard is a quantitative grade standard of 20-HETE-d, (item No. 390030) that has been prepared specifically for mass spectrometry and related applications where quantitative reproducibility is required. The solution has been prepared gravimetrically and is supplied in a deactivated glass ampule sealed under argon. The concentration was verified by comparison to an independently prepared calibration standard. This 20-HETE-d₆ MaxSpec[®] standard is guaranteed to meet identity, purity, stability, and concentration specifications and is provided with a batch-specific certificate of analysis. Ongoing stability testing is performed to ensure the concentration remains accurate throughout the shelf life of the product. Note: The amount of solution added to the vial is in excess of the listed amount. Therefore, it is necessary to accurately measure volumes for preparation of calibration standards. Follow recommended storage and handling conditions to maintain product quality.

References

- 1. McGiff, J.C. and Quilley, J. 20-HETE and the kidney: Resolution of old problems and new beginnings. Am. J. Physiol. 277(3 Pt 2), R607-R623 (1999).
- 2. Gebremedhin, D., Lange, A.R., Lowry, T.F., et al. Production of 20-HETE and its role in autoregulation of cerebral blood flow. Circ. Res. 87(1), 60-65 (2000).
- 3 Prakash, C., Zhang, J.Y., Falck, J.R., et al. 20-Hydroxyeicosatetraenoic acid is excreted as a glucuronide conjugate in human urine. Biochem. Biophys. Res. Commun. 185(2), 728-733 (1992).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 02/13/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM