

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



MEGA-10

Item No. 25703

CAS Registry No.:	85261-20-7	
Formal Name:	1-deoxy-1-[methyl(1-oxodecyl) amino]-D-glucitol	
Synonym:	N-decanoyl-N-Methylglucamine	ОН ОН
MF:	C ₁₇ H ₃₅ NO ₆	
FW:	349.5	HO, \downarrow \uparrow \checkmark \downarrow \checkmark \checkmark \checkmark \checkmark \checkmark
Purity:	≥98%	он он о
Supplied as:	A crystalline solid	
Storage:	-20°C	
Stability:	≥2 years	
Information represents	s the product specifications. Batch specifi	c analytical results are provided on each certificate of analysis.

Laboratory Procedures

MEGA-10 is supplied as a crystalline solid. A stock solution may be made by dissolving the MEGA-10 in the solvent of choice. MEGA-10 is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF), which should be purged with an inert gas. The solubility of MEGA-10 in ethanol is approximately 1 mg/ml and approximately 30 mg/ml in DMSO and DMF.

MEGA-10 is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, MEGA-10 should first be dissolved in DMSO and then diluted with the aqueous buffer of choice. MEGA-10 has a solubility of approximately 0.33 mg/ml in a 1:2 solution of DMSO:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

MEGA-10 is a nonionic detergent that can be used to solubilize membrane proteins.¹ It has a critical micelle concentration (CMC) of 4.88 mM under no salt conditions, which decreases in the presence of sodium chloride.² MEGA-10 has been used to reconstitute amino acid transporters from rat liver plasma membrane vesicles into artificial phospholipid membranes.³ It has also been used to solubilize the melibiose transport carrier from E. coli membranes and reconstitute it into liposomes.⁴

References

- 1. Hildreth, J.E.K. N-D-Gluco-N-methylalkanamide compounds, a new class of non-ionic detergents for membrane biochemistry. Biochem. J. 207(2), 363-366 (1982).
- 2. Molina-Bolívar, J.A., Hierrezuelo, J.M., and Carnero Ruiz, C. Self-assembly, hydration, and structures in N-decanoyl-N-methylglucamide aqueous solutions: Effect of salt addition and temperature. J. Colloid. Interface Sci. 313(2), 656-664 (2007).
- 3. Quesada, A.R. and McGivan, J.D. A rapid method for the functional reconstitution of amino acid transport systems from rat liver plasma membranes. Partial purification of System A. Biochem J. 255(3), 963-969 (1988).
- 4. Hanatani, M., Nishifuji, K., Futai, M., et al. Solubilization and reconstitution of membrane proteins of Escherichia coli using alkanoyl-N-methylglucamides. J. Biochem. 95(5), 1349-1353 (1984).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 10/11/2018

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM