

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



1-Tridecanoyl-rac-glycerol

Item No. 26961

CAS Registry No.:	71958-74-2	
Formal Name:	tridecanoic acid,	
Synonyms:	2,3-dihydroxypropyl ester MG(13:0/0:0/0:0), 1-Monotridecanoin	
MF:	C ₁₆ H ₃₂ O ₄	
FW:	288.4	
Purity:	≥95%	
Supplied as:	A solid	
Storage:	-20°C	
Stability:	≥2 years	
Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.		

Laboratory Procedures

1-Tridecanoyl-rac-glycerol is supplied as a solid. A stock solution may be made by dissolving the 1-tridecanoyl-rac-glycerol in the solvent of choice, which should be purged with an inert gas. 1-Tridecanoylrac-glycerol is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of 1-tridecanoyl-rac-glycerol in these solvents is approximately 15, 20, and 25 mg/ml, respectively.

1-Tridecanoyl-rac-glycerol is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, 1-tridecanoyl-rac-glycerol should first be dissolved in DMF and then diluted with the aqueous buffer of choice. 1-Tridecanoyl-rac-glycerol has a solubility of approximately 0.3 mg/ml in a 1:2 solution of DMF:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

1-Tridecanoyl-rac-glycerol is a monoacylglycerol that contains tridecanoic acid (Item No. 19725) at the sn-1 position. It has been used as an internal standard for the quantification of mono- and diacylglycerols in buttermilk.1

Reference

1. Fagan, P., Wijesundera, C., and Watkins, P. Determination of mono- and di-acylglycerols in milk lipids. J. Chromatogr. A. 1054(1-2), 251-259 (2004).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 03/13/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM