

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

# SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# PRODUCT INFORMATION



## CTOP (trifluoroacetate salt)

Item No. 27377

Formal Name: D-phenylalanyl-L-cysteinyl-L-tyrosyl-

> D-tryptophyl-L-ornithyl-L-threonyl-3mercapto-L-valyl-L-threoninamide-cyclic (2→7)-disulfide, trifluoroacetate salt

Synonym: Phe-Cys-Tyr-Trp-Orn-Thr-Pen-Thr MF: C<sub>50</sub>H<sub>67</sub>N<sub>11</sub>O<sub>11</sub>S<sub>2</sub>• XCF<sub>3</sub>COOH

FW: 1,062.3 Purity:

λ<sub>max</sub>: 220, 269, 386 nm UV/Vis.:

Supplied as: A solid -20°C Storage: Stability: ≥2 years

XCF<sub>2</sub>COOH

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

#### **Laboratory Procedures**

CTOP (trifluoroacetate salt) is supplied as a solid. A stock solution may be made by dissolving the CTOP (trifluoroacetate salt) in the solvent of choice, which should be purged with an inert gas. CTOP (trifluoroacetate salt) is slightly soluble in ethanol.

CTOP (trifluoroacetate salt) is slightly soluble in aqueous solutions. To enhance aqueous solubility, dilute the organic solvent solution into aqueous buffers or isotonic saline. If performing biological experiments, ensure the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. We do not recommend storing the aqueous solution for more than one day.

### Description

CTOP is a peptide antagonist of the  $\mu$ -opioid receptor (IC<sub>50</sub> = 2.8 nM).<sup>1</sup> It is selective for  $\mu$ -opioid receptors over  $\delta$ -opioid receptors (IC<sub>50</sub> = 13,500 nM) and somatostatin receptors (IC<sub>50</sub> = 24,000 nM). CTOP inhibits the antinociceptive effect of morphine in the tail flick test in mice ( $ED_{50} = 0.018$  nmol) and reverses morphine-induced increases in locomotor activity (ED $_{50}$  = 0.02 nmol).<sup>2</sup> It also induces withdrawal symptoms in a mouse model of chronic, but not acute, morphine dependence.

#### References

- 1. Pelton, J.T., Kazmierski, W., Gluya, K., et al. Design and synthesis of conformationally constrained somatostatin analogues with high potency and specificity for μ opioid receptors. J. Med. Chem. 29(11), 2370-2375 (2013).
- 2. Gulya, K., Kriván, M., Nyolczas, N., et al. Central effects of the potent and highly selective μ opioid antagonist D-Phe-Cys-Tyr-D-Trp-Orn-Thr-Pen-Thr-NH2 (CTOP) in mice. Eur. J. Pharmacol. 150(3), 355-360 (1988).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

## WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 07/24/2019

### **CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM