



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

PRODUCT INFORMATION



Long-chain Fatty Acid Methyl Ester Mixture (porcine ovary)

Item No. 28662

Supplied as: A liquid

Storage: -20°C

Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.


Description and Data

Long-chain fatty acid methyl ester mixture contains 14 omega fatty acids and three saturated fatty acids isolated from porcine ovaries and supplied as methyl esters. It is intended for use as a reference standard for the detection of these long-chain fatty acids from various sample types.

Included are:

| Elution Order* | Item Name |
|----------------|---|
| 1 | Myristic Acid methyl ester |
| 2 | Palmitic Acid methyl ester |
| 3 | Palmitoleic Acid methyl ester |
| 4 | Stearic Acid methyl ester |
| 5 | Oleic Acid methyl ester |
| 6 | Vaccenic Acid methyl ester |
| 7 | Linoleic Acid methyl ester |
| 8 | Jacaric Acid methyl ester |
| 9 | Linolenic Acid methyl ester |
| 10 | 11(Z)-Eicosenoic Acid methyl ester |
| 11 | 11(Z),14(Z)-Eicosadienoic Acid methyl ester |
| 12 | Dihomo-γ-Linolenic Acid methyl ester |
| 13 | Arachidonic Acid methyl ester |
| 14 | Eicosapentaenoic Acid methyl ester |
| 15 | 7(Z),10(Z),13(Z),16(Z)-Docosatetraenoic Acid methyl ester |
| 16 | Docosapentaenoic Acid methyl ester |
| 17 | Docosahexaenoic Acid methyl ester |

*GC Column Used: Omegawax 250, 30 m x 0.25 mm x 0.25 μm



WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM

WWW.CAYMANCHEM.COM