

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



.NO₂

Dantrolene-d₄ Item No. 29451

Formal Name: 1-[[[5-(4-nitrophenyl)-furanyl]

methylenelaminol-2,4-

imidazolidinedione

MF: $C_{14}H_6D_4N_4O_5$

FW: 318.3

Chemical Purity: ≥98% (Dantrolene)

Deuterium

 \geq 99% deuterated forms (d₁-d₄); \leq 1% d₀ Incorporation:

A solid Supplied as: -20°C Storage: Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Dantrolene- d_{Δ} is intended for use as an internal standard for the quantification of dantrolene (Item No. 14326) by GC- or LC-MS. The accuracy of the sample weight in this vial is between 5% over and 2% under the amount shown on the vial. If better precision is required, the deuterated standard should be quantitated against a more precisely weighed unlabeled standard by constructing a standard curve of peak intensity ratios (deuterated versus unlabeled).

Dantrolene-d₄ is supplied as a solid. A stock solution may be made by dissolving the dantrolene-d₄ in the solvent of choice, which should be purged with an inert gas. Dantrolene-d₁ is soluble in DMSO.

Description

Dantrolene is an inhibitor of calcium release from the sarcoplasmic reticulum ($IC_{50} = 0.3 \mu M$).¹ It binds to sarcoplasmic reticulum vesicles isolated from normal and malignant hyperthermia-susceptible (MHS) pigs with K_d values of 0.3 and 0.09 μ M, respectively. Dantrolene reduces spontaneous calcium wave frequency and amplitude in the presence of calmodulin in isolated mouse cardiomyocytes with IC₅₀ values of 0.42 and 0.19 μM, respectively.² It reduces the magnitude of electrically stimulated twitch tensions in isolated rat extensor digitorum longus and soleus muscles (IC₅₀ = 3 μg/ml) and MHS pigs (ED₅₀ = 0.85 mg/kg).^{3,4} Formulations containing dantrolene have been used in the treatment of malignant hyperthermia.

References

- 1. Kobayashi, S., Yano, M., Suetomi, T., et al. Dantrolene, a therapeutic agent for malignant hyperthermia, markedly improves the function of failing cardiomyocytes by stabilizing interdomain interactions within the ryanodine receptor. J. Am. Coll. Cardiol. 53(21), 1993-2005 (2009).
- Y.W., Gomez-Hurtado, N., Walweel, K., et al. Essential role of calmodulin in RyR inhibition by dantrolene. Mol. Pharmacol. 88(1), 57-63 (2015).
- Kotsias, B.A. and Muchnik, S. Reversible effect of dantrolene sodium on twitch tension of rat skeletal muscle. Arch. Neurol. 35(4), 234-236 (1978).
- Nelson, T.E. and Flewellen, E.H. Rationale for dantrolene vs. procainamide for treatment of malignant hyperthermia. Anesthesiology 50(2), 118-122 (1979).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 01/21/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM