

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# PRODUCT INFORMATION



СООН

# Palmitic Acid-do MaxSpec® Standard

Item No. 30355

CAS Registry No.: 1173022-49-5 Formal Name: hexadecanoic-

13,13,14,14,15,15,16,16,16-do acid

Synonyms: C16:0-d<sub>9</sub>. Cetylic Acid-d<sub>9</sub>

 $C_{16}H_{23}D_{9}O_{2}$ MF: FW: 265.5 **Purity:** ≥95%

Supplied as: A solution in chloroform; in a deactivated glass ampule

Concentration: 100 μg/ml (nominal); see certificate of analysis for verified concentration

Storage:

Stability: ≥1 year; Stability testing is ongoing to ensure concentration accuracy. The certificate of analysis and

product expiry date will be updated upon completion of testing.

Special Conditions: Store upright and unopened at -20°C. Warm to room temperature prior to opening.

Light sensitive.

## Description

Palmitic acid-do is intended for use as an internal standard for the quantification of palmitic acid (Item No. 10006627) by GC- or LC-MS. Palmitic acid is a common 16-carbon saturated fat that represents 10-20% of human dietary fat intake and comprises approximately 25 and 65% of human total plasma lipids and saturated fatty acids, respectively. $^{1,2}$  Acylation of palmitic acid to proteins facilitates anchoring of membrane-bound proteins to the lipid bilayer and trafficking of intracellular proteins, promotes protein-vesicle interactions, and regulates various G protein-coupled receptor functions. Red blood cell palmitic acid levels are increased in patients with metabolic syndrome compared to patients without metabolic syndrome and are also increased in the plasma of patients with type 2 diabetes compared to individuals without diabetes.<sup>3,4</sup>

Palmitic acid-d<sub>o</sub> MaxSpec<sup>®</sup> standard is a quantitative grade standard of palmitic acid (Item No. 30380) that has been prepared specifically for mass spectrometry or any application where quantitative reproducibility is required. The solution has been prepared gravimetrically and is supplied in a deactivated glass ampule sealed under argon. The concentration was verified by comparison to an independently prepared calibration standard. The verified concentration is provided on the certificate of analysis. This palmitic acid-do MaxSpec® standard is guaranteed to meet identity, purity, stability, and concentration specifications and is provided with a batch-specific certificate of analysis. Ongoing stability testing is performed to ensure the concentration remains accurate throughout the shelf life of the product. Note: The amount of solution added to the vial is in excess of the listed amount. Therefore, it is necessary to accurately measure volumes for preparation of calibration standards. Follow recommended storage and handling conditions to maintain product quality.

## References

- 1. Fatima, S., Hu, X., Gong, R.-H., et al. Palmitic acid is an intracellular signaling molecule involved in disease development. Cell. Mol. Life Sci. 76(13), 2547-2557 (2019).
- Santos, M.J., López-Jurado, M., Llopis, J., et al. Influence of dietary supplementation with fish oil on plasma fatty acid composition in coronary heart disease patients. Ann. Nutr. Metab. 39(1), 52-62 (1995).
- Yi, L.-Z., He, J., Liang, Y.-Z., et al. Plasma fatty acid metabolic profiling and biomarkers of type 2 diabetes mellitus based on GC/MS and PLS-LDA. FEBS Lett. 580(30), 6837-6845 (2006).
- Kabagambe, E.K., Tsai, M.Y., Hopkins, P.N., et al. Erythrocyte fatty acid composition and the metabolic syndrome: A National Heart, Lung, and Blood Institute GOLDN study. Clin. Chem. 54(1), 154-162 (2008).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information Buyer agrees to purchase the mater can be found on our website.

Copyright Cayman Chemical Company, 06/02/2021

## **CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM