

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Nordeoxycholic Acid

Item No. 30837

CAS Registry No.: 53608-86-9

Formal Name: (3α,5β,12α)-3,12-dihydroxy-24-

norcholan-23-oic acid

Synonyms: 23-NDCA, 23-Nordeoxycholic Acid

MF: $C_{23}H_{38}O_4$ FW: 378.6 **Purity:** ≥98% Supplied as: A solid Storage: -20°C Stability: ≥2 years

Ĥ HO

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Nordeoxycholic acid is supplied as a solid. A stock solution may be made by dissolving the nordeoxycholic acid in the solvent of choice, which should be purged with an inert gas. Nordeoxycholic acid is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of nordeoxycholic acid in ethanol and DMSO is approximately 10 mg/ml and approximately 15 mg/ml in DMF.

Nordeoxycholic acid is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, nordeoxycholic acid should first be dissolved in DMF and then diluted with the aqueous buffer of choice. Nordeoxycholic acid has a solubility of approximately 0.2 mg/ml in a 1:4 solution of DMF:PBS (pH 7.2) using this method. We do not recommend storing the aqueous solution for more than one day.

Description

Nordeoxycholic acid is a metabolite of the bile acid norcholic acid and a 23-carbon derivative of deoxycholic acid (Item Nos. 20756 | 18231). Levels of nordeoxycholic acid are decreased in the liver of rats in a high-fat diet model of non-alcoholic fatty liver disease (NAFLD) and serum levels are lower in men compared with women.^{2,3} Nordeoxycholic acid has commonly been used as an internal standard for the quantification of bile acids in various sample types by GC- and LC-MS.^{4,5}

References

- 1. Nakatomi, F., Kihira, K., Kuramoto, T., et al. Intestinal absorption and metabolism of norcholic acid in rats. J. Pharmacobiodyn 8(7), 557-563 (1985).
- 2. Tang, Y., Zhang, J., Li, J., et al. Turnover of bile acids in liver, serum and caecal content by high-fat diet feeding affects hepatic steatosis in rats. Biochim. Biophys. Acta Mol. Cell Biol. Lipids 1864(10), 1293-1304
- 3. Xie, G., Wang, Y., Wang, X., et al. Profiling of serum bile acids in a healthy Chinese population using UPLC-MS/MS. J. Proteome Res. 14(2), 850-859 (2015).
- 4. Bergeron, A., Furtado, M., and Garofolo, F. Importance of using highly pure internal standards for successful liquid chromatography/tandem mass spectrometric bioanalytical assays. Rapid Commun. Mass Spectrom. 23(9), 1287-197 (2009).
- 5. Evrard, E. and Janssen, G. Gas-liquid chromatographic determination of human fecal bile acids. J. Lipid Res. **9(2)**, 226-236 (1968).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 07/21/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM