

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATIO



Glycocholic Acid MaxSpec® Standard

Item No. 31351

CAS Registry No.: 475-31-0

Formal Name: N- $[(3\alpha,5\beta,7\alpha,12\alpha)-3,7,12-trihydroxy-24-$

oxocholan-24-yl]-glycine

Synonyms: Cholylglycine, GCA

MF: $C_{26}H_{43}NO_{6}$ FW: 465.6 **Purity:** >95%

Supplied as: A solution in methanol; in a deactivated glass ampule

Concentration: 100 μg/ml (nominal); see certificate of analysis for verified concentration

Storage:

Stability: ≥2 years; Stability testing is ongoing to ensure concentration accuracy. The certificate of analysis and

product expiry date will be updated upon completion of testing.

Special Conditions: Store upright and unopened at -20°C. Warm to room temperature prior to opening.

Light sensitive.

Description

Glycocholic acid is a glycine-conjugated form of the primary bile acid cholic acid (Item No. 20250) and has roles in the emulsification of fats. 1.2 It reduces expression of the gene encoding the farnesoid X receptor (FXR) and increases expression of the genes encoding the bile acid receptors TGR5 and S1PR2 in SNU-245 cells when used at a concentration of 1.6 μmol/ml.³ Glycocholic acid (250 μM) increases the intracellular accumulation and cytotoxicity of epirubicin (Item No. 12091) in Caco-2 cells, as well as decreases expression of the genes encoding multidrug resistance protein 1 (MDR1), MDR-associated protein 1 (MRP1), and MRP2 when used alone or in combination with epirubicin.⁴ It increases absorption of epirubicin into everted sacs of rat ileum and jejunum when used at a concentration of 250 µM. The bile acid composition ratio of glycocholic acid is elevated in the bile of patients with cholangiocarcinoma compared with patients with pancreatic cancer or benign biliary diseases.³ Serum levels of glycocholic acid are elevated in patients with hepatocellular carcinoma compared with healthy individuals.²

Glycocholic acid MaxSpec® standard is a quantitative grade standard of glycocholic acid (Item No. 20276) that has been prepared specifically for mass spectrometry or any application where quantitative reproducibility is required. The solution has been prepared gravimetrically and is supplied in a deactivated glass ampule sealed under argon. The concentration was verified by comparison to an independently prepared calibration standard. This glycocholic acid MaxSpec® standard is guaranteed to meet identity, purity, stability, and concentration specifications and is provided with a batch-specific certificate of analysis. Ongoing stability testing is performed to ensure the concentration remains accurate throughout the shelf life of the product. **Note:** The amount of solution added to the vial is in excess of the listed amount. Therefore, it is necessary to accurately measure volumes for preparation of calibration standards. Follow recommended storage and handling conditions to maintain product quality.

References

- 1. Lefebyre, P., Cariou, B., Lien, F., et al. Physiol. Rev. 89(1), 147-191 (2009).
- Guo, C., Xie, C., Ding, P., et al. J. Chromatogr. B. Analyt. Technol. Biomed. Life Sci. 1072, 315-319 (2018).
- 3. Song, W.-S., Park, H.-M., Ha, J.M., et al. Sci. Rep. 8(1), 11088 (2018).
- 4. Lo, Y.L., Ho, C.T., and Tsai, F.L. Eur. J. Pharm. Sci. 35(1-2), 52-67 (2008).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information Buyer agrees to purchase the mater can be found on our website.

Copyright Cayman Chemical Company, 05/20/2021

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA **PHONE:** [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM