

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Taurocholic Acid-d₄ MaxSpec[®] Standard

Item No. 31375

CAS Registry No.: 252030-90-3

Formal Name: $2-[[(3\alpha,5\beta,7\alpha,12\alpha)-3,7,12-trihydroxy-$

24-oxocholan-24-yl-2,2,4,4-d₄]amino]-

ethanesulfonic acid

Synonyms: Ethanesulfonic Acid-d₄, TCA-d₄

MF: $C_{26}H_{41}D_4NO_7S$

FW: 518.7 **Purity:** ≥95%

Supplied as: A solution in methanol; in a deactivated glass ampule

Concentration: 100 μg/ml (nominal); see certificate of analysis for verified concentration

Storage: -20°C

Stability: ≥2 years; Stability testing is ongoing to ensure concentration accuracy. The certificate of analysis and

product expiry date will be updated upon completion of testing.

Special Conditions: Store upright and unopened at -20°C. Warm to room temperature prior to opening.

Light sensitive.

Description

Taurocholic acid-d₄ (TCA-d₄) is intended for use as an internal standard for the quantification of TCA (Item No. 16215) by GC- or LC-MS. TCA is a taurine-conjugated form of the primary bile acid cholic acid (Item No. 20250). Serum levels of TCA are decreased in patients with Crohn's disease and those with ulcerative colitis with no extraintestinal manifestations but are increased in patients with ulcerative colitis accompanied by hepatobiliary disease.²

TCA-d₄ MaxSpec® standard is a quantitative grade standard of TCA-d₄ (Item No. 21891) that has been prepared specifically for mass spectrometry or any application where quantitative reproducibility is required. The solution has been prepared gravimetrically and is supplied in a deactivated glass ampule sealed under argon. The concentration was verified by comparison to an independently prepared calibration standard. The verified concentration is provided on the certificate of analysis. This TCA-d₄ MaxSpec[®] standard is guaranteed to meet identity, purity, stability, and concentration specifications and is provided with a batch-specific certificate of analysis. Ongoing stability testing is performed to ensure the concentration remains accurate throughout the shelf life of the product. Note: The amount of solution added to the vial is in excess of the listed amount. Therefore, it is necessary to accurately measure volumes for preparation of calibration standards. Follow recommended storage and handling conditions to maintain product quality.

References

- 1. Lefebvre, P., Cariou, B., Lien, F., et al. Role of bile acids and bile acid receptors in metabolic regulation. Physiol. Rev. 89(1), 147-191 (2009).
- 2. Gnewuch, C., Liebisch, G., Langmann, T., et al. Serum bile acid profiling reflects enterohepatic detoxification state and intestinal barrier function in inflammatory bowel disease. World J. Gastroenterol. **15(25)**, 3134-3141 (2009).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information Buyer agrees to purchase the m can be found on our website.

Copyright Cayman Chemical Company, 01/11/2022

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM