

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Oxotremorine (sesquifumarate)

Item No. 35707

CAS Registry No.: 17360-35-9

1-[4-(1-pyrrolidinyl)-2-butyn-1-yl]-2-Formal Name:

pyrrolidinone, 2E-butenedioate (2:3)

MF: C₁₂H₁₈N₂O • 3/2C₄H₄O₄

FW: 380.4 **Purity:** ≥95% Supplied as: A solid Storage: -20°C Stability: ≥4 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Oxotremorine (sesquifumarate) is supplied as a solid. A stock solution may be made by dissolving the oxotremorine (sesquifumarate) in the solvent of choice, which should be purged with an inert gas. Oxotremorine (sesquifumarate) is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide. The solubility of oxotremorine (sesquifumarate) in these solvents is approximately 10, 14, and 16 mg/ml, respectively.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of oxotremorine (sesquifumarate) can be prepared by directly dissolving the solid in aqueous buffers. The solubility of oxotremorine (sesquifumarate) in PBS (pH 7.2) is approximately 3 mg/ml. We do not recommend storing the aqueous solution for more than one day.

Description

Oxotremorine is an agonist of muscarinic acetylcholine receptors (mAChRs).¹ It binds to M₁, M₂, M₃, and M_A mAChRs (K_is = 3.31, 9.33, 5.24, and 7.08 μ M, respectively, for the human receptors). It inhibits acetylcholine (ACh; İtem No. 23829) release in rat cerebral cortical slices when used at a concentration of 10 μM, an effect that can be reversed by the non-selective mAChR antagonist atropine (Item No. 12008).² Oxotremorine induces tremors, salivation, and lacrimation in mice (ED_{50} s = 0.11, 0.22, and 0.28 mg/kg).³ lt induces analgesia in the hot plate test in mice when administered at a dose of 0.5 mg/kg.4 Oxotremorine has been used to induce parkinsonian symptoms as a model of Parkinson's disease in mice.

References

- 1. Jakubík, I., Bacáková, L., El-Fakahany, E.E., et al. Positive cooperativity of acetylcholine and other agonists with allosteric ligands on muscarinic acetylcholine receptors. Mol. Pharmacol. 52(1), 172-179 (1997).
- Szerb, J.C. and Somogyi, G.T. Depression of acetylcholine release from cerebral cortical slices by cholinesterase inhibition and by oxotremorine. Nat. New Biol. 241(108), 121-122 (1973).
- Moon, M.W., Chidester, C.G., Heier, R.F., et al. Cholinergic activity of acetylenic imidazoles and related compounds. J. Med. Chem. 34(8), 2314-2327 (1991).
- Pleuvry, B.J. and Tobias, M.A. Comparison of the antinociceptive activities of physostigmine, oxotremorine and morphine in the mouse. Br. J. Pharmacol. 43(4), 706-714 (1971).
- Castensson, S., Sievertsson, H., Lindeke, B., et al. Studies on the inhibition of oxotremorine induced tremor by a melanocyte-stimulating hormone release-inhibiting factor, thyrotropin releasing hormone and related peptides. FEBS Lett. 44(1), 101-105 (1974).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 12/12/2022

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM