

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Thymine

Item No. 38749

CAS Registry No.: 65-71-4

Formal Name: 5-methyl-2,4(1H,3H)-pyrimidinedione Synonyms: 5-Methyluracil, NSC 14705, NSC 168663

MF: $C_5H_6N_2O_2$ FW: 126.1 **Purity:** ≥98%

UV/Vis.: λ_{max} : 210, 264 nm

A solid Supplied as: -20°C Storage: Stability: ≥4 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



Laboratory Procedures

Thymine is supplied as a solid. A stock solution may be made by dissolving the thymine in the solvent of choice, which should be purged with an inert gas. Thymine is soluble in organic solvents such as ethanol, DMSO, and dimethyl formamide (DMF). The solubility of thymine in ethanol is approximately 2 mg/ml and approximately 20 mg/ml in DMSO and DMF.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. Organic solvent-free aqueous solutions of thymine can be prepared by directly dissolving the solid in aqueous buffers. Thymine is slightly soluble in PBS (pH 7.2). We do not recommend storing the aqueous solution for more than one day.

Description

Thymine is a pyrimidine base. 1 It forms complementary base pairs with the purine adenine (Item No. 18148) in DNA.¹ Thymine is produced by the catabolism of thymidine (Item No. 20519) via thymidine phosphorylase.² It is replaced with uracil (Item No. 26088) in RNA.³

References

- 1. Ghannam, J.Y., Wang, J., and Jan, A. Biochemistry, DNA Structure. StatPearls Publishing, 1-5 (2021).
- Brown, N.S. and Bicknell, R. Thymidine phosphorylase, 2-deoxy-D-ribose and angiogenesis. Biochem. J. 334(Pt 1), 1-8 (1998).
- 3. Banoub, J.H., Newton, R.P., Esmans, E., et al. Recent developments in mass spectrometry for the characterization of nucleosides, nucleotides, oligonucleotides, and nucleic acids. Chem. Rev. 105(5), 1869-1915 (2005).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 07/26/2023

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM