



**SZABO  
SCANDIC**

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](http://linkedin.com/company/szaboscandic)



# PRODUCT INFORMATION



## OP-145 (trifluoroacetate salt)

Item No. 38841

Formal Name:	N-acetyl-L-isoleucylglycyl-L-lysyl-L- $\alpha$ -glutamyl-L-phenylalanyl-L-lysyl-L-arginyl-L-isoleucyl-L-valyl-L- $\alpha$ -glutamyl-L-arginyl-L-isoleucyl-L-lysyl-L-arginyl-L-phenylalanyl-L-leucyl-L-arginyl-L- $\alpha$ -glutamyl-L-leucyl-L-valyl-L-arginyl-L-prolyl-L-leucyl-L-argininamide, trifluoroacetate salt	Ac—Ile—Gly—Lys—Glu—Phe—Lys—Arg—Ile—Val—Glu—Arg—Ile—Lys—Arg—Phe—Leu—Arg—Glu—Leu—Val—Arg—Pro—Leu—Arg—NH <sub>2</sub>
Synonyms:	acetyl-IGKEFKRIVERIKRFLRELVRPLR-amide, P60.4-Ac	• XCF <sub>3</sub> COOH
Peptide Sequence:	Ac-IGKEFKRIVERIKRFLRELVRPLR-NH <sub>2</sub>	
MF:	C <sub>142</sub> H <sub>246</sub> N <sub>46</sub> O <sub>31</sub> • XCF <sub>3</sub> COOH	
FW:	3,093.8	
Purity:	≥95%	
Supplied as:	A solid	
Storage:	-20°C	
Stability:	≥4 years	

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

### Laboratory Procedures

OP-145 (trifluoroacetate salt) is supplied as a solid. A stock solution may be made by dissolving the OP-145 (trifluoroacetate salt) in water. We do not recommend storing the aqueous solution for more than one day.

### Description

OP-145 is a synthetic antimicrobial peptide and derivative of the antimicrobial peptide LL-37 (Item No. 24461).<sup>1</sup> It is active against *E. coli*, *P. aeruginosa*, *C. albicans*, and *A. niger* (MICs = 2, 3, 6, and 18  $\mu$ M, respectively). Poly(lactic-co-glycolic acid) (PLGA) microspheres containing OP-145 prevent biofilm formation in a mouse model of bone infection induced by surgical femoral fracture followed by bacterial inoculation at the fracture site.<sup>2</sup>

### References

1. Nell, M.J., Tjabringa, G.S., Wafelman, A.R., et al. Development of novel LL-37 derived antimicrobial peptides with LPS and LTA neutralizing and antimicrobial activities for therapeutic application. *Peptides* **27**(4), 649-660 (2006).
2. Cheng, Y., Qin, J., Huang, Y., et al. The antimicrobial effects of PLGA microspheres containing the antimicrobial peptide OP-145 on clinically isolated pathogens in bone infections. *Sci. Rep.* **12**(1), 14541 (2022).

**WARNING**  
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

**SAFETY DATA**

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

**WARRANTY AND LIMITATION OF REMEDY**

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 07/31/2023

**CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD  
ANN ARBOR, MI 48108 - USA

**PHONE:** [800] 364-9897  
[734] 971-3335

**FAX:** [734] 971-3640

CUSTSERV@CAYMANCHEM.COM  
WWW.CAYMANCHEM.COM