

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



Product Information



Arachidonic Acid-d_o

Item No. 390010

CAS Registry No.: 69254-37-1

Formal Name: 5Z,8Z,11Z,14Z-eicosatetraenoic-

5,6,8,9,11,12,14,15-d_g acid

Synonym:

MF: $C_{20}H_{24}D_8O_2$ FW: 312.5

Chemical Purity:

≥96% Arachidonic Acid

Deuterium

Incorporation: \geq 99% deuterated forms (d₁-d₈); \leq 1% d₀

Stability: ≥1 year at -20°C

Supplied as: A solution in methyl acetate

Laboratory Procedures

Arachidonic acid-d₈ contains eight deuterium atoms at the 5, 6, 8, 9, 11, 12, 14, and 15 positions. It is intended for use as an internal standard for the quantification of arachidonic acid by GC- or LC-mass spectrometry (MS). For long term storage, we suggest that arachidonic acid-d₈ be stored as supplied at -20°C. It should be stable for at least one year.

Arachidonic acid-d₈ is supplied as a solution in methyl acetate. To change the solvent, simply evaporate the methyl acetate under a gentle stream of nitrogen and immediately add the solvent of choice. Solvents such as ethanol, DMSO, or dimethyl formamide purged with an inert gas can be used. The solubility of arachidonic acid-d₈ in these solvents is approximately 100 mg/ml.

Arachidonic acid- d_8 is used as an internal standard for the quantification of arachidonic acid by stable isotope dilution MS. The accuracy of the sample weight in this vial is between 5% over and 2% under the weight indicated on the vial. If better precision is required, the deuterated standard should be quantitated against a more precisely weighed unlabeled standard arachidonic acid by constructing a standard curve of peak intensity ratios (deuterated versus unlabeled).

Arachidonic acid is an essential fatty acid and a precursor for all prostaglandins, thromboxanes, and leukotrienes. Virtually all cellular arachidonic acid is esterified in membrane phospholipids where its presence is tightly regulated through multiple interconnected pathways. Free arachidonic acid is a transient, critical substrate for the biosynthesis of eicosanoid second messengers. Receptor-stimulated release, metabolism, and re-uptake of free arachidonate are all important aspects of cell signaling and inflammation.²

References

- 1. Nixon, A.B., Greene, D.G., and Wykle, R.L. Comparison of acceptor and donor substrates in the CoA-independent transacylase reaction in human neutrophils. Biochim. Biophys. Acta 1300, 187-196 (1996).
- Burgoyne, R.D. and Morgan, A. The control of free arachidonic acid levels. Trends Biochem. Sci. 15, 365-366 (1990).

Related Products

For a list of related products please visit: www.caymanchem.com/catalog/390010

WARNING: This product is for laboratory research only: not for administration to humans. Not for human or veterinary DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Cayman Chemical Company makes **no warranty or guarantee** of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman **warrants only** to the original customer that the material will <u>meet our specifications</u>

purpose, sustainity and inestinancement, which care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, includental or consequential damages, even if Cayman is informed about their possible existence.

This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees.

Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a <u>refund</u> of the purchase price, or at Cayman's option, the <u>replacement</u>, at no cost to Buyer, of all material that

Buyer's exclusive remedy and Caymans sole liability nereunder snail be limited to a ferund of the purchase price, of at Cayman within constitutions.

Said refund or replacement is conditioned on Buyer giving written notice to Cayman within thirty (30) days after arrival of the material at its destination. Failure of Buyer to give said notice within thirty (30) days shall constitute a waiver by Buyer of all claims hereunder with respect to said material.

For further details, please refer to our Warranty and Limitation of Remedy located on our website and in our catalog.

Copyright Cayman Chemical Company, 03/04/2015

Cayman Chemical

Mailing address

1180 E. Ellsworth Road Ann Arbor, MI 48108 USA

Phone

(800) 364-9897 (734) 971-3335

(734) 971-3640

custserv@caymanchem.com

www.caymanchem.com