

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

PRODUCT INFORMATION



C2 Ceramide (d18:1/2:0)

Item No. 62510

| CAS Registry No.: Formal Name: | 3102-57-6 N-[(1S,2R,3E)-2-hydroxy-1- (hydroxymethyl)-3-heptadecen-1- yl]-acetamide | | H N |
|-----------------------------------|--|-------------------------------------|--------------|
| Synonyms: | C2 Ceramide, Ceramide (d18:1/2:0), Cer(18:1/2:0), N-acetoyl-D- <i>erythro</i> -sphingosine | | ОН |
| MF: | C ₂₀ H ₃₉ NO ₃ | | он |
| FW: | 341.5 | | |
| Purity: | ≥98% | | |
| Supplied as: | A crystalline solid | | |
| Storage: | -20°C | | |
| Stability: | As supplied, 1 year from the QC data stored properly | e provided on the Certificate of An | alysis, when |

Laboratory Procedures

Concentrated stock solutions of C2 ceramide (d18:1/2:0) can be prepared by dissolving the crystalline solid in an organic solvent such as ethanol, methanol, chloroform, or DMSO. The solubility of C2 ceramide (d18:1/2:0) in these solvents is approximately 10 mg/ml. Stock solutions will be stable for six months stored at -20°C. Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Also, ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. We do not recommend storing the aqueous solution for more than one day.

Description

C2 ceramide is a biologically active, cell permeable, and less hydrophobic analog of natural ceramides. It has been shown to induce cell differentiation and inhibit cell growth in HL-60 cells.^{2,3} C2 ceramide activates a cytosolic protein phosphatase in T9 cells at levels as low as 0.1 μ M.⁴

References

- 1. Obeid, L.M., Linardic, C.M., Karolak, L.A., et al. Programmed cell death induced by ceramide. Science 259, 1769-1771 (1993).
- 2. Okazaki, T., Bielawska, A., Bell, R.M., et al. Role of ceramide as a lipid mediator of 1α , 25-dihydroxyvitamin D₃-induced HL-60 cell differentiation. J. Biol. Chem. 265, 15823-15831 (1990).
- 3. Bielawska, A., Linardic, C.M., and Hannun, Y.A. Ceramide-mediated biology. Determination of structural and stereospecific requirements through the use of N-acyl-phenylaminoalcohol analogs. J. Biol. Chem. 267, 18493-18497 (1992).
- 4. Dobrowsky, R.T. and Hannun, Y.A. Ceramide stimulates a cytosolic protein phosphatase. J. Biol. Chem. 267, 5048-5051 (1992).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 02/19/2018

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM