

PRODUCT INFORMATION



COX-FA Positive Control

Item No. 700110

Overview

Contents: This vial contains 50 µl of ovine COX-1.
Storage: -80°C
Stability: As supplied, 2 years from the QC date provided on the Certificate of Analysis, when stored properly

Laboratory Procedures

COX-FA Positive Control should be kept on ice when thawed. Prior to performing the COX Fluorescent Activity Assay Kit (Item No. 700200), dilute 10 µl of enzyme with 490 µl of diluted Assay Buffer (Item No. 760154) and store on ice. The diluted enzyme is stable for one hour.

Precaution

This product may not perform as described if used with other assay reagents or protocols.

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the [complete](#) Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 08/16/2023

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897
[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM