



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

PRODUCT INFORMATION



FFA Assay Cofactor Mixture 1

Item No. 700312

Overview

Contents: The vial contains a lyophilized mixture of acyl coenzyme A synthetase and ascorbate oxidase.

Stability: ≥ 1 year at -20°C

Laboratory Procedures

Prior to use in the Free Fatty Acid Fluorometric Assay Kit (Item No. 700310), reconstitute one vial with 6 ml of FFA Assay Buffer (Item No. 700311). Place vial on ice. The reconstituted mixture is stable for three hours on ice. The mixture will be combined with the reconstituted FFA Assay Cofactor Mixture 2 (Item No. 700304) prior to utilizing in the assay. One vial is enough to assay 60 wells. If additional wells are being utilized, then reconstitute a second vial. *NOTE: Ascorbate oxidase eliminates ascorbic acid interference.*

Cautions

For laboratory research use only. This reagent is not for human or veterinary use. This product has been tested and formulated to work exclusively with the Free Fatty Acid Fluorometric Assay Kit from Cayman (Item No. 700310). This product may not perform as described if used with other assay reagents or protocols.

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent *via* email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 09/08/2015

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897
[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM