

# Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

# **Product Information**



### SOD Assay Sample Buffer (10X)

Item No. 706003

**Contents:** This vial contains 500 mM Tris-HCl, pH 8.0. Stability: ≥1 year at -20°C

#### Laboratory Procedures

Prior to use in the Superoxide Dismutase Assay Kit (Item No. 706002), dilute 2 ml of Sample Buffer concentrate with 18 ml of HPLC-grade water for assaying 96 wells. Prepare additional Sample Buffer as needed. This final Sample Buffer (50 mM Tris-HCl, pH 8.0) should be used to prepare the SOD Standard (Item No. 706005) and dilute the SOD Assay Xanthine Oxidase (Item No. 706006) and SOD samples prior to assaying in the Superoxide Dismutase Assay Kit (Item No. 706002). When stored at 4°C, this diluted Sample Buffer is stable for at least six months.

#### CAUTIONS

For laboratory research use only. This reagent is not for human or veterinary use. This product has been tested and formulated to work exclusively with the Superoxide Dismutase Assay Kit from Cayman (Item No. 706002). This product may not perform as described if used with other assay reagents or protocols.

#### **Related Products**

DIAGNOSTIC OR THERAPEUTIC USE.

WARRANTY AND LIMITATION OF REMEDY

SAFETY DATA

at the time of delivery.

For a list of related products please visit: www.caymanchem.com/catalog/706003

#### **Cayman Chemical**

#### **Mailing address**

1180 E. Ellsworth Road Ann Arbor, MI 48108 USA

Phone (800) 364-9897 (734) 971-3335

Fax (734) 971-3640

#### E-Mail

custserv@caymanchem.com

#### Web www.caymanchem.com

ar the time of delivery. Cayman will carry out its delivery obligations with due care and skill. Thus, in no event will Cayman have **any obligation or liability**, whether in tort (including negligence) or in contract, for any direct, indirect, indirect incidental or consequential damages, even if Cayman is informed about their possible existence. This limitation of liability does not apply in the case of intentional acts or negligence of Cayman, its directors or its employees. Buyer's **exclusive remedy** and Cayman's sole liability hereunder shall be limited to a <u>refund</u> of the purchase price, or at Cayman's option, the <u>replacement</u>, at no cost to Buyer, of all material that

This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, or inhale. Do not get in eyes, on skin, or on clothing. Wash thoroughly after handling. This information contains some, <u>but not all</u>, of the information required for the safe and proper use of this material. Before use, the user <u>must</u> review the <u>complete</u> Safety Data Sheet, which has been sent *via* email to your institution.

Cayman Chemical Company makes no warranty or guarantee of any kind, whether written or oral, expressed or implied, including without limitation, any warranty of fitness for a particular purpose, suitability and merchantability, which extends beyond the description of the chemicals hereof. Cayman warrants only to the original customer that the material will meet our specifications

WARNING: This product is for laboratory research only: not for administration to humans. Not for human or veterinary

Buyer's exclusive remedy and Laymans sole hability neterinder shall be infinited to a terminal of the particular process of an experimental strength of the particular process of the particular proces of the