



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 


PRODUCT INFORMATION



Nervonoyl Ethanolamide

Item No. 9001746

CAS Registry No.: 887405-21-2
Formal Name: N-(2-hydroxyethyl)-15Z-tetracosenamide
MF: $C_{26}H_{51}NO_2$
FW: 409.7
Purity: $\geq 98\%$
Supplied as: A crystalline solid
Storage: $-20^{\circ}C$
Stability: ≥ 4 years



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Nervonoyl ethanolamide is supplied as a crystalline solid. A stock solution may be made by dissolving the nervonoyl ethanolamide in the solvent of choice, which should be purged with an inert gas. Nervonoyl ethanolamide is soluble in the organic solvent ethanol at a concentration of approximately 5 mg/ml.

Description

Nervonoyl ethanolamide is a member of the family of fatty N-acyl ethanolamines collectively called endocannabinoids.¹⁻⁴ The relative importance of this ethanolamine metabolite has not been determined.

References

1. Bachur, N.R. and Udenfriend, S. Microsomal synthesis of fatty acid amides. *J. Biol. Chem.* **241**(6), 1308-1313 (1966).
2. Doetsch, P.W., Zastawny, T.H., Martin, A.M., *et al.* Monomeric base damage products from adenine, guanine, and thymine induced by exposure of DNA to ultraviolet radiation. *Biochemistry* **34**(3), 737-742 (1995).
3. Saghatelian, A., Trauger, S.A., Want, E.J., *et al.* Assignment of endogenous substrates to enzymes by global metabolite profiling. *Biochemistry* **43**(45), 14332-14339 (2004).
4. Buczynski, M.W., Svensson, C.I., Dumlao, D.S., *et al.* Inflammatory hyperalgesia induces essential bioactive lipid production in the spinal cord. *J. Neurochem.* **114**(4), 981-993 (2010).

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the [complete](#) Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 04/30/2024

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897
[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM