



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## NORMAL GOAT SERUM

**CL1200**  
**CL1200-10**  
**CL1200-100**  
**CL1200-500**

### DESCRIPTION:

Cedarlane normal goat serum is derived from a closed herd of adult goats. It has been processed for minimum haemolysis and filtered to 0.22µm. CL1200 is preserved with 0.1% sodium azide to allow for long term storage at 4°C.

### PRESENTATION:

CL1200: 10 ml, liquid with 0.1% NaN<sub>3</sub>.  
CL1200-10: 10 ml, frozen, filtered 0.22 µm, sterile.  
CL1200-100: 100ml, frozen, filtered 0.22 µm, sterile.  
CL1200-500: 500ml, frozen, filtered 0.22 µm, sterile.

### STORAGE/STABILITY:

Store CL1200 at 4°C.

**Note: Sera may precipitate in storage.**

Store CL1200-10, CL1200-100 and CL1200-500 at -20°C or below.

### APPLICATIONS:

Blocking sera, control sera, a serum additive in dilution buffers.

**FOR RESEARCH USE ONLY**

rev/da/apr/00

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)