



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.
Please contact CEDARLANE® for lot specific information.

GOAT BLOOD IN ELSEVER'S SOLUTION

CL1200-100A
CL1200-500A
CL1200-1000A

DESCRIPTION:

A fifty percent (v/v) solution of horse blood in Alsever's solution is obtained aseptically, close to the time of shipment

PRESENTATION:

100 ml, liquid
500 ml, liquid
1000ml, liquid

STORAGE/STABILITY:

Store at +4°C. Expiry is three weeks from the date of collection, which is indicated on the label.

Cedarlane employs highest quality standards in the collection of goat blood.

Our goat blood is guaranteed sterile, and is tested for sterility.

FOR RESEARCH USE ONLY

EJ/04/05

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com