



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

NORMAL SHEEP SERUM

CL1300-50

CL1300-100

CL1300-500

DESCRIPTION:

Cedarlane Laboratories has developed a normal sheep serum specifically for use in immunohistopathology. It has been processed for minimum haemolysis and filtered to 0.22µm.

PRESENTATION:

CL1300-50: 50 ml, frozen, filtered 0.22µm, sterile.

CL1300-100: 100 ml, frozen, filtered 0.22µm, sterile.

CL1300-500: 500 ml, frozen, filtered 0.22µm, sterile.

STORAGE/STABILITY:

Store CL1300 at -20°C or below.

APPLICATIONS:

Blocking sera, control sera, a serum additive in dilution buffers.

FOR RESEARCH USE ONLY

rev/lh/nov/98

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com