



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

BOVINE BLOOD IN ELSEVER'S SOLUTION

CL1700-100A

CL1700-500A

CL1700-1000A

DESCRIPTION:

A fifty percent (v/v) solution of bovine blood in Alsever's solution is obtained aseptically, close to the time of shipment.

PRESENTATION:

100 ml, liquid

500 ml, liquid

1000 ml, liquid

STORAGE AND STABILITY:

Store at 4°C. Expiry is three weeks from the date of collection that is indicated on the label.

Cedarlane employs highest quality standards in the collection of bovine blood.

Our bovine blood is guaranteed sterile, and is tested for sterility.

FOR RESEARCH USE ONLY

rev/lh/nov/98

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com