



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

*Please contact CEDARLANE® for lot specific information.*

## Sheep Anti-Rabbit Thrombin Purified IgG

**CL20121A-RB**

**Lot #:**

**Exp:**

**Immunogen:** Rabbit  $\alpha$ -Thrombin inactivated with PPACK.

**Presentation:** **1.0 ml** of immunoglobulin purified from sheep antiserum. The IgG concentration is **10 mg/ml** in a solution of 10mM HEPES, pH 7.4, 0.15 M NaCl and 50% (v/v) glycerol.

**Storage:** Store between -10 and -20°C. Mix gently before use.

**Specificity:** IgG prepared from antiserum was analysed by immunoelectrophoresis and found to react with purified thrombin and prothrombin. Only one arc was seen with normal plasma, no arc was seen with barium-adsorbed rabbit plasma.

**Application:** This product can be used in immunoassays for prothrombin and thrombin as well as thrombin-inhibitor complexes.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

**For Research Use Only**

JV1/19/01

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



**toll free: 1-800-268-5058**  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)