



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

*Please contact CEDARLANE® for lot specific information.*

## SWINE ANTI-RABBIT IMMUNOGLOBULIN FLORESCEIN CONJUGATED

**CL6200F**

**LOT: 435**

### **DESCRIPTION:**

The antiserum is prepared by immunizing swine with rabbit IgG. The antiserum is specific for rabbit immunoglobulin: It does not react with rabbit serum proteins other than immunoglobulins: it does not cross-react with immunoglobulins of other species.

### **PRESENTATION:**

FITC - conjugated antibody in 0.05M carbonate buffer.

**FORMAT:** 1 ml lyophilized.

### **STORAGE AND RECONSTITUTION:**

Store at -20°C or below before reconstitution. Reconstitute with 1 ml of cold distilled water. Aliquot and store at -70° after reconstitution. Prepare aliquots in appropriate volume for single use to avoid repeated freeze-thawing.

### **LOT SPECIFICATIONS:**

Lot: 435

Protein Concentration: 8.0 ± 0.2 mg / ml

F:P Ratio 8.8

Recommended working dilution 1:20-1:80

**FOR RESEARCH USE ONLY**

SV/02/22/99

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



*toll free: 1-800-268-5058*  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)