



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

*Please contact CEDARLANE® for lot specific information.*

## NORMAL RAT SERUM

**CL7000**

**CL7000-50**

**CL7000-500S**

### DESCRIPTION:

Cedarlane Laboratories has developed a normal rat serum specifically for use in immunohistopathology. It has been processed for minimum haemolysis and filtered to 0.45µm non-sterile, 0.22µm for sterile products.

### PRESENTATION:

CL7000: 5 ml, frozen, filtered 0.45µm, non-sterile.

CL7000-50: 50ml, frozen, filtered 0.45µm, non-sterile.

CL7000-500S: 500ml, frozen, filtered 0.22µm, sterile.

### STORAGE/STABILITY:

Store at -20°C or below.

### APPLICATIONS:

Blocking sera, control sera, a serum additive in dilution buffers.

**FOR RESEARCH USE ONLY**

rev/DA/Oct/00

For more information or to place an order please contact...

**CEDARLANE®**  
**LABORATORIES LIMITED**



*toll free: 1-800-268-5058*  
*in North America*

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: [www.cedarlanelabs.com](http://www.cedarlanelabs.com) • e-mail: [info@cedarlanelabs.com](mailto:info@cedarlanelabs.com)