



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

**FITC Anti-Mouse Macrophage (F4/80)
Monoclonal Antibody**

CL8940F
CL8940F-3
LOT: 14011206

DESCRIPTION:

Cedarlane's anti-mouse F4/80 monoclonal antibody reacts with the mouse macrophage F4/80 antigen, which is a 160 kD plasma membrane component on mouse mononuclear phagocytes. The F4/80 antigen is found on most macrophages, and on macrophage precursors from M-CFC onward. Expression of this antigen is increased upon maturation. F4/80 is found in low levels on activated macrophages and eosinophils. Dendritic leukocytes may be negative or express F4/80 in low levels.

Applications include flow cytometry^{7,8,9,11}. This clone is also reported to work in immunohistochemistry, both frozen and paraffin sections^{5,12}, and ELISA¹¹.

PRESENTATION:

100 µg (CL8940F) or 300 µg (CL8940F-3) FITC conjugated Ig buffered in PBS, 0.1% sodium azide (NaN₃) and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/ml.

STORAGE/STABILITY:

Store at 4°C. For long term storage, aliquot and freeze unused portion at -20°C in volumes appropriate for single usage. Avoid freeze/thaw cycles and prolonged exposure to light. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

Continued Overleaf..

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

SPECIFICATIONS:

Clone: C1:A3-1

Specificity: Mouse Macrophage (F4/80)

Ig Class: Rat IgG_{2b}

Antibody Concentration: 0.1 mg/ml

FLOW CYTOMETRY ANALYSIS:

Method:

1. Prepare a cell suspension in media A. For cell preparations, deplete the red blood cell population with Lympholyte[®]-M cell separation medium (CL5030).
2. Wash 2 times.
3. Resuspend the cells to a concentration of 2×10^7 cells/ml in media A. Add 50 μ l of this suspension to each tube (each tube will then contain 1×10^6 cells, representing 1 test).
4. To each tube, add 1.0 μ g* of **CL8940F or CL8940F-3** per 10^6 cells.
5. Vortex the tubes to ensure thorough mixing of antibody and cells.
6. Incubate the tubes for 30 minutes at 4°C.
(It is recommended that the tubes are protected from light, since most fluorochromes are light sensitive.)
7. Wash 2 times at 4°C.
8. Resuspend the cell pellet in 50 μ l ice cold media B.
9. Transfer to suitable tubes for flow cytometric analysis containing 15 μ l of propidium iodide at 0.5 mg/ml in PBS. This stains dead cells by intercalating in DNA.

Media:

- A. Phosphate buffered saline (pH 7.2) + 5% normal serum of host species + sodium azide (100 μ l of 2M sodium azide in 100 mls).
- B. Phosphate buffered saline (pH 7.2) + 0.5% Bovine serum albumin + sodium azide (100 μ l of 2M sodium azide in 100 mls).

Results:

Tissue Distribution by Flow Cytometry Analysis:
(Representative Histogram)


Mouse Strain: BALB/c

Cell Concentration : 1×10^6 cells per test

Antibody Concentration Used: 1.0 μ g/ 10^6 cells

Isotypic Control: FITC Rat IgG_{2b}

Continued Overleaf..



Cell Source: Peritoneal Macrophage

Percentage of cells stained above control: 84.0%

N.B. Appropriate control samples should always be included in any labeling studies.

*** For optimal results in various applications, it is recommended that each investigator determine dilutions appropriate for individual use.**

REFERENCES:

1. Austyn, J. M., Gordon, S. (1981) F4/80, a monoclonal antibody directed specifically against the mouse macrophage. *Eur. J. Immuno.* 11: 805-815.
2. Berard, J., *et al.* (1997). Abnormal regulation of retinoic acid receptor γ_2 Expression and compromised allograft rejection in transgenic mice expressing antisense sequences to retinoic acid receptor γ_1 and γ_3 . *Jour. of Immuno.* 159: 2586-2598.
3. Gordon, S. *et al.* (1992). Macrophages in tissue and in vitro. *Current opinion in immunology.* 4: 25-32.
4. Gordon, S. *et al.* (1992). Antigen markers of macrophage differentiation in murine tissues. *Current topics in microbio. and immuno.* 181: 1-37.
5. Handbook of Experimental Immunology, Ed. Weir, D.M. (Chapter 43).
6. Hume, D.A., Perry, V.H., Gordon, S. "The Mononuclear Phagocyte System of the Mouse Defined by Immunohistochemical Localizations of Antigen F4/80: Macrophages Associated with Epithelia", *The Anatomical Record* (1984), Volume 210, pp 503-512.
7. Leenen, P., *et al.* (1998). Heterogeneity of mouse spleen Dendritic cells: in vivo phagocytic activity, expression of macrophage markers, and subpopulation turnover. *Jour. of Immuno.* 160: 2166-2173.
8. Leenen, P. *et al.* (1986). Murine macrophage cell lines can be ordered in a linear differentiation sequence. *Differentiation.* 32:157-164.
9. Li, X. Q., *et al.* (1998). Immunohistochemical detection of testicular macrophage during the period of postnatal maturation in the mouse. *Int. J. Androl.* 21(6): 370-376.

Continued Overleaf...

10. McKnight, A. J., *et al.* (1996). Molecular cloning of F4/80, a murine macrophage-restricted cell surface glycoprotein with homology to the G-protein-linked transmembrane 7 hormone receptor family. *J. Biol. Chem.* 271(1): 486-489.
11. Szu-Hee Lee, Starky, P.M., Gordon, S. "Quantitative Analysis of Total Macrophage Content in Adult Mouse Tissues", *J. Exp. Med.* (1985), Volume 161 pp 475-489.
12. Whiteland, J.L *et al* (1995). Immunohistochemical detection of T cell subsets and other leukocytes in paraffin embedded rat and mouse tissues with monoclonal antibodies .*J. Histochem. Cytochem.* 43: 313-320.

FOR RESEARCH USE ONLY

® is a Registered Trademark of Cedarlane Laboratories Limited.

CH 03/03/05