



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Rabbit Anti-Ochratoxin A, Total IgG

CLAS11-1683

- Background:** Ochratoxin A, a toxin produced by *Aspergillus ochraceus* and *Penicillium verrucosum*, is one of the most abundant food-contaminating mycotoxins in the world. The toxin has been found in the tissues and organs of animals, including human blood and breast milk. Ochratoxin A can cause immunosuppression and immunotoxicity in animals.
- Format:** Rabbit polyclonal, total IgG, lyophilized in PBS pH 7.4.
- Immunogen:** BSA-conjugated ochratoxin A.
- Quantity:** 2.5 mg
- Reconstitution:** For reconstitution add 2.5 ml of sterile water.
- Storage:** Store lyophilized and reconstituted product at 4°C.
- Application:** ELISA (ELISA)
- Recommended dilution:** The optimal working dilution should be determined by the investigator.
- Confirmed Reactivity:** *Aspergillus ochraceus*
- Not reactive in:** No confirmed exceptions from predicted reactivity are currently known.
- Additional information:** Purified antibodies were dialyzed against deionized water.
- Cross reaction with ochratoxin B - 2%.

JK 01/25/22

**Laboratory Reagent For Research Use Only**

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2015 and ISO 13485:2016  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)