



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Turkey Red Blood Cell Suspension (10%)

CLC1180

Lot: 16118004

DESCRIPTION:

The whole blood is aseptically collected and prepared to order from our colony of animals. These washed cell suspensions are useful for the titration of complement, absorption procedures, testing for agglutinins, and for the preparation of stroma as particulate reagents.

SIZE:

50 mL, liquid

PRESENTATION:

Supplied as a washed and freshly prepared 10% suspension of Turkey Red Blood Cells in Alsever's solution.

STORAGE:

Store at 4°C. **DO NOT FREEZE.**

STABILITY:

Short shelf life. For optimal viability, the blood cells should be used within 7 days of shipment date.

APPLICATIONS:

This reagent is suitable for use in absorption, agglutination and hemolytic procedures.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

Laboratory Reagent For Research Use Only

XT 02/18/16

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com