



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Sheep Red Blood Cell Suspension (10%)

**CLC6800**  
**Lot: 16036811**

### **DESCRIPTION:**

The whole blood is aseptically collected and prepared to order from fresh sheep blood in Alsevers. These washed cell suspensions are useful for the titration of complement, absorption procedures, testing for agglutinins, and for the preparation of stroma as particulate reagents.

### **SIZE:**

100 mL, liquid

### **PRESENTATION:**

Supplied as a washed and freshly prepared 10% suspension of Sheep Red Blood Cells in Alsever's solution.

### **STORAGE:**

Store at 4°C. **DO NOT FREEZE.**

### **STABILITY:**

Short shelf life. For optimal viability, the blood cells should be used within 7 days of shipment date.

### **APPLICATIONS:**

This reagent is suitable for use in absorption, agglutination and hemolytic procedures.

**Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory**

**Laboratory Reagent For Research Use Only**

XT 10/11/15

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)