



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

## Purified Goat IgG Isotype Control

**CLCC10200**

**DESCRIPTION:** Purified Goat IgG

**VOLUME:** 2.0 ml

**ANTIBODY CONTENT:** 2.0 mg

**PREPARATION:**

Normal Goat IgG fractionated from goat serum.

**APPLICATIONS:**

This isotype control is provided in liquid format for application in flow cytometry. 1 µg is the maximum amount of reagent required for  $1 \times 10^6$  cells. Lesser amount of reagent may be required and it is recommended that the customer determine the optimum amount for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

**BUFFERS AND PRESERVATIVES:**

This product is provided in phosphate buffered saline and it contains no preservatives.

**STORAGE/STABILITY:**

Stable at 4°C. For long term storage, aliquot and freeze unused portions at -20°C in volumes appropriate for single usage. Avoid repeated freeze/thaw cycles. If the reagent is being diluted, it is recommended that the quantity to be used within one week be diluted.

**Laboratory Reagent For Research Use Only**

AF 11/01/10

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

**In CANADA: Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

**In the USA: Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)