



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

PE Goat IgG Isotype Control

CLCC10204

LOT:

DESCRIPTION:

Normal Goat IgG fractionated from goat serum was conjugated with R-Phycoerythrin.

VOLUME: 1.0 ml

ANTIBODY CONTENT: 100 tests

PREPARATION:

Antibodies were affinity purified by salt-based precipitation of the immunoglobulin from normal goat serum. The purity of the IgG fraction is verified by ELISA using a panel of anti-goat antibodies. The resulting antibody was then conjugated to PE to provide optimal fluorochrome to protein ratio for flow cytometry. Any free dye was removed by column chromatography.

APPLICATIONS:

This isotype control is provided in liquid format for application in flow cytometry. 10 µl is the maximum amount of reagent required for 1×10^6 cells. Lesser amount of reagent may be required and it is recommended that the customer determine the optimum amount for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

BUFFERS AND PRESERVATIVES:

Product contains 0.1% sodium azide. A highly purified grade of BSA has been added as a stabilizing protein to bring the final concentration to 4-5mg/ml after conjugation.

STORAGE/STABILITY:

Stable at 4°C. **DO NOT FREEZE.** Light exposure should be avoided. If the reagent is being diluted, it is recommended that the quantity to be used within one week be diluted. Check label for expiry date.

Laboratory Reagent For Research Use Only

BA 05/18/12

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com