



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®]
www.cedarlanelabs.com



Delivering You Today's Innovations for
the Science of Tomorrow™

**Purified Goat F(ab')₂ IgG
Isotype Control**

CLCC11300

LOT:

DESCRIPTION: Purified Goat F(ab')₂ isotype control.

PREPARATION & APPLICATION:

Normal goat IgG was digested with pepsin and any aggregates, Fc or undigested goat IgG were removed by column chromatography. This isotype control is provided in liquid format for application in flow cytometry. 1 µg is the maximum amount of reagent required for 1 x 10⁶ cells. Lesser amount of reagent may be required and it is recommended that the customer determine the optimum amount for each application.

VOLUME: 1.0 ml

ANTIBODY CONTENT: 5 mg

BUFFERS AND PRESERVATIVES:

This product is presented in 100 mM borate buffered saline, pH 8.2. No preservatives have been added.

STORAGE/STABILITY:

Stable at 4°C. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

Laboratory Reagent For Research Use Only

JV 06/24/09

Visit our website for your local distributor.

CEDARLANE[®]



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com