



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

FITC Mouse IgG_{2b} Isotypic Control

CLCMG2B01

LOT:

DESCRIPTION: FITC Mouse IgG_{2b} isotypic control

SIZE: 0.5 ml

ANTIBODY CONTENT: 50 µg

FLUOROCHROME/PROTEIN RATIO: 4-7

APPLICATION: Each lot is tested in flow cytometry using human peripheral blood lymphocytes (PBL). From this testing it is recommended that 10 µl of antibody be used per 1 x 10⁶ cells. Because results may vary, lesser amounts of reagent may be sufficient and it is recommended that the customer determine the optimum amount for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

BUFFERS AND PRESERVATIVES: This reagent is provided in phosphate buffered saline (PBS) and contains 0.1% sodium azide (NaN₃) as a preservative. A highly purified grade of BSA has been added as a stabilizing agent to bring to final protein concentration to 4-5 mg/ml after conjugation.

Continued Overleaf..

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information
website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

STORAGE/STABILITY: Stable at 4°C. For long term storage, aliquot and freeze unused portions at -20°C in volumes appropriate for single use. Avoid repeated freeze/thaw cycles and prolonged exposure to light. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

FOR RESEARCH USE ONLY

CH 04/18/05