



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®]
www.cedarlanelabs.com



Delivering You Today's Innovations for
the Science of Tomorrow™

**PE Mouse IgG_{2b}
Isotype Control**

CLCMG2B04
LOT: 2B510807

DESCRIPTION: PE Mouse IgG_{2b} isotype control

CLONE: A-1

VOLUME: 0.5 ml

ANTIBODY CONTENT: 50 µg

APPLICATIONS:

For use in flow cytometry, 10 µl is the maximum amount of reagent required for 1x10⁶ cells. Lesser amounts of reagent may be sufficient and it is recommended that the customer determine the optimum amount of product for each application. Each lot is tested in flow cytometry using freshly harvested mouse splenocytes and/or human peripheral blood lymphocytes (PBL). When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

BUFFERS & PRESERVATIVES:

This reagent is provided in phosphate buffered saline (PBS) and contains 0.1% sodium azide (NaN₃) as a preservative, and a stabilizing agent.

STORAGE / STABILITY:

Store at 4°C. **DO NOT FREEZE.** If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

FOR RESEARCH USE ONLY

AF 07/27/09

Visit our website for your local distributor.

CEDARLANE[®]



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com