



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

PE-Cy5 Mouse IgG_{2b} Isotype Control

CLCMG2B06

LOT:

DESCRIPTION: PE-Cy5 Mouse IgG_{2b} isotype control

VOLUME: 0.5 ml

ANTIBODY CONTENT: 50 µg

APPLICATIONS: For use in flow cytometry, 10 µl is the maximum amount of reagent required for 10⁶ cells. Lesser amounts of reagent may be sufficient and it is recommended that the customer determine the optimum amount of product for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

PE-Cy5: This conjugate can be used as a third colour reagent in flow cytometry. It requires a 650 nm long pass filter in the FL-3 channel. FL2 - FL3 compensation will be in the range of 1%.

BUFFERS & PRESERVATIVES: This reagent is provided in phosphate buffered saline (PBS) and it contains 0.1% sodium azide (NaN₃) as a preservative. A highly purified grade of BSA has been added as a stabilizing protein to bring the final concentration to 4-5 mg/ml after conjugation.

Continued Overleaf..

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com

STORAGE / STABILITY: Store at 4°C. **DO NOT FREEZE.** Check label for expiry date. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

FOR RESEARCH USE ONLY

CH 04/18/05