



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

TECHNICALLY *Speaking*

Place your order with CEDARLANE® or your local distributor.

Please contact CEDARLANE® for lot specific information.

Rabbit anti-GFAP

CLCN0021

Lot:

Product:

Rabbit anti-GFAP, affinity purified pAbs against glial fibrillary acid protein (GFAP), an astrocyte marker of central nervous system. Human, mouse and rat specific.

Size:

50 µg in 50 µl PBS

Immunogen:

GFAP fusion protein

Suggested Usage:

The product can be used at 1:1,000 to 1:2,000 for immunoblot testing or 1:500 – 1:1000 for IHC or ICC. The product can be used in combination with Cedarlane's conjugates such as CLDRAP, CLDRCY3, CLDRFITC or CLGRHRP.

Short Term Storage:

2-8°C for one week

Long Term Storage:

-20°C as an undiluted liquid for up to two years.

Note:

For in vitro research use only, not for diagnostic or therapeutic use.

TN 09/12/05

For more information or to place an order please contact...

CEDARLANE®
LABORATORIES LIMITED



toll free: 1-800-268-5058
in North America

phone: (905) 878-8891 • fax: (905) 878-7800

5516 - 8th Line, R.R.#2, Hornby, Ontario, CANADA L0P 1E0

or visit our website for a list of our international distributors including contact information

website: www.cedarlanelabs.com • e-mail: info@cedarlanelabs.com