



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

**PE Rat IgG_{1κ}
Isotypic Control**

CLCR104
LOT: 13041213

DESCRIPTION: PE Rat IgG_{1κ} isotypic control

CLONE: KLH/G1-2-2

SIZE: 1.0 ml

ANTIBODY CONTENT: 100 µg

APPLICATION:

This isotype control is provided in liquid form for use in flow cytometry. 0.2 µg is the maximum amount of reagent required for 1×10^6 cells. Lesser amounts of reagent may be sufficient and it is recommended that the customer determine the optimum amount for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

BUFFERS & PRESERVATIVE:

This reagent is provided in phosphate buffered saline and contains 0.02% sodium azide (NaN₃) as a preservative, and a stabilizing agent.

STORAGE/STABILITY:

Stable at 4°C. **DO NOT FREEZE.** Avoid prolonged exposure to light. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

Laboratory Reagent For Research Use Only

BA 12/27/13

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com