



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

FITC Rat IgG_{2b} Isotypic Control

CLCR2B01

DESCRIPTION: FITC Rat IgG_{2b} isotypic control

CLONE: RTK4530

SIZE: 0.5 ml

ANTIBODY CONTENT: 50 µg

APPLICATION:

This isotype control is provided in liquid form for use in flow cytometry. 1 µg is the maximum amount of reagent required for 1x10⁶ cells. Lesser amounts of reagent may be sufficient, and it is recommended that the customer determine the optimum amount for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

BUFFERS & PRESERVATIVE:

This reagent is provided in phosphate buffered saline and contains 0.09% sodium azide (NaN₃) as a preservative. A highly purified grade of BSA has been added as a stabilizing protein to bring the final protein concentration to 4-5 mg/ml.

STORAGE/STABILITY:

Stable at 4°C. For long term storage, aliquot and freeze unused portions at -20°C in volumes appropriate for single use. Avoid repeated freeze/thaw cycles and prolonged exposure to light. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted.

Laboratory Reagent For Research Use Only

MW 01/20/22

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com