



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically  
Speaking

**CEDARLANE**<sup>®</sup>   
[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

Conveniently Delivering You Today's Innovations  
for the Science of Tomorrow™

**APC Rat IgG<sub>2b</sub>  
Isotypic Control**

**CLCR2B05**  
**LOT:** 13051211

**DESCRIPTION:** APC Rat IgG<sub>2b</sub> isotypic control

**SIZE:** 1.0 ml

**PROTEIN CONTENT:** 100 µg

**APPLICATION :**

This isotype control is provided in liquid form for use in flow cytometry. 10 µl is the maximum amount of reagent required for 1x10<sup>6</sup> cells. Lesser amounts of reagent may be sufficient, and it is recommended that the customer determine the optimum amount for each application. When using isotype controls as negative controls, it is recommended that the amount of isotype control and the testing antibody used be equivalent.

**BUFFERS & PRESERVATIVE:**

This reagent is presented in phosphate buffered saline containing 0.1% sodium azide (NaN<sub>3</sub>) as preservative. A highly purified grade of BSA has been added as a stabilizing protein to bring the final concentration to 4-5 mg/ml after conjugation.

**STORAGE/STABILITY:**

Stable at 4°C. **DO NOT FREEZE.** Avoid prolonged exposure to light. If the reagent is being diluted, it is recommended that only the quantity to be used within one week be diluted. Check label for expiry date.

**Laboratory Reagent For Research Use Only**

BA 12/19/11

Visit our website for your local distributor.

**CEDARLANE**<sup>®</sup> 

[www.cedarlanelabs.com](http://www.cedarlanelabs.com)

An ISO 9001:2000 and ISO 13485:2003  
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020  
e-mail: [general@cedarlanelabs.com](mailto:general@cedarlanelabs.com)

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138  
e-mail: [service@cedarlanelabs.com](mailto:service@cedarlanelabs.com)