



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Recombinant Human Heparanase-1 Active

CLENZ1032
CLENZ1032-2
CLENZ1032-3

Introduction: Heparanase is an endo β -D-glucuronidase, which degrades heparan sulfate side chains of heparan sulfate proteoglycans (HSPGs) in the extracellular matrix. Heparanase plays an important role in ECM degradation, facilitating the migration and extravasation of tumor cells and inflammatory leukocytes (1,2,3). Upon degradation, heparanase releases growth factors and cytokines that stimulate cell proliferation and chemotaxis (4,5). Heparanase is a heterodimer comprised of a 50 kDa subunit harboring the active site and a 8 kDa subunit. It is produced as a latent 65 kDa precursor and proteolytically processed to its active form (1,6). Heparanase is highly expressed in myeloid leukocytes (i.e. neutrophils) in platelets and in human placenta. Human heparanase was found to be upregulated in various types of primary tumors, correlating in some cases with increased tumor invasiveness and vascularity and with poor prospective survival (7,8).

Description: Heparanase Active Enzyme is produced in CHO cells. The protein is purified by several orthogonal chromatography steps.

Source: CHO cells.

Presentation: 1 μ g (CLENZ1032), 4 μ g (CLENZ1032-2), or 16 μ g (CLENZ1032-3). Heparanase Active Enzyme is supplied in 20mM Acetate buffer and 750mM NaCl pH 5.4.

Stability: Store at -80°C , avoid repeated freeze-thaw cycles.

Purity: Greater than 95.0% as determined by:

- (a) Analysis by RP-HPLC.
- (b) Analysis by SDS-PAGE.

Specificity: Heparanase Active Enzyme is identified by Western blot analysis with polyclonal rabbit anti-HPA1 antibodies as 2 subunits of 8-kDa and 50-kDa.

Biological Activity: The specific activity of Heparanase Active Enzyme in-house standard is about 0.7 Units (1 unit = 1 μ mole of reducing ends of heparan sulfate substrate formed per minute per mg Heparanase Active Enzyme at 37°C). The enzymatic activity of each Heparanase Active Enzyme batch is comparable to the standard as determined by activity assay in which immobilized heparan, released due to heparanase activity, is quantified colorimetrically. Recommended reaction buffer: 20 mM Citrate Phosphate buffer, pH 5.4; 50mM NaCl; 1mM CaCl_2 .

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com

References:

1. Vlodaysky, I. et al., Nat. Med. 5 (1999) 793–802.
2. Vlodaysky, I and Friedman, Y., J. Clin. Invest. 108 (2001) 341–347.
3. Parish, C.R. et al., Biochem. Biophys. Acta 1471 (2001) M99–M108.
4. Vlodaysky, I. et al., Cancer Metastasis Rev. 9, (1990) 203-226.
5. Bashkin, P. et al., Biochemistry 28 (1989), 1737-1743.
6. Fairbanks, M.B. et al., J. Biol. Chem. 274 (1999) 29587– 29590.
7. Koliopanos, A. et al., Cancer Res. 61 (1992) 4655–4659.
8. Gohji, K. et al., Int. J. Cancer 95 (2001) 295–301.

Laboratory Reagent For Research Use Only

JK 09/07/21